



UNIVERSITAT DE
BARCELONA

**PLASTIC SCINTILLATORS
and
RELATED ANALYTICAL PROPOSALS
for
RADIONUCLIDE ANALYSIS**

José F. García, Alex Tarancón, Héctor Bagán.

Departament d'Enginyeria Química i Química Analítica

LSC2017 Conference, 1-5 May, Copenhagen



Begining of Plastic Scintillators

The Theory and Practice of Scintillation Counting

J.B. Birks

Pergamon Press 1964

Liquid Scintillation and Plastic Scintillation starts 70's

Liquid Scintillation

Focus: radionuclide analysis

→ LS Counter – LS Spectrometers

→ Scintillation Cocktails

→ LSC success: . Sample preparation for measurements (easy and homogeneous)

. Many applications

Plastic Scintillation

Focus: large area detectors, dosimetry

Radionuclide analysis

Plastic Scintillators and related analytical proposals

- Plastic Scintillators composition and sample preparation
- Scintillation mechanism
- Direct radionuclide determinations
- Selective radionuclide determinations
- Future challenges


Plastic Scintillators: composition and sample preparation

Composition

- . Solvent: . Polystyrene, Polyvinyltoluene → linear (water insoluble)
 - . Divinylbenzene → crosslinked (aggressive and organic insoluble)
- . Secondary solvent: Naphtalene, Disopropilnaphtalene
- . Scintillators: PPO, POPOP, p-T, bis-MSB

Format

- . Plastic scintillator microspheres (PSm) (10 – 300 μm)
- . Plastic scintillator foils (PSf) (50 – 100 μm)
- . Plastic scintillator pellets (3 mm)



Plastic Scintillators: composition and sample preparation

Sample preparation for measurement.

- Steps:
- . Plastic scintillator
 - . Sample solution
 - . Homogenization


LS vs PS: . PS time ↑

- . PS difficulty to homogenize
- . PS expensive than LS cocktails
- . PS sample solution and PSm can be segregated (no mixed waste)
- . PS sample stability (no phase separation)




Scintillation mechanism

Scintillation mechanism




Quenching

LS



Scintillation mechanism


Scintillation mechanism



Quenching




LS

PS




Scintillation mechanism

Scintillation



Quenching



PS


Particle Quenching $\rightarrow \phi \uparrow$ EFF \downarrow
 Optical Quenching $\rightarrow \phi \uparrow$ EFF \uparrow

PSm ϕ 10 – 300 μ m


Direct radionuclide determinations - PSm

Alpha emitter radionuclides.

	EFF	
PSm	^{241}Am (%) (5552 keV)	^{236}U (%) (4573 keV)
D (70 μm)	96	75
D300 (300 μm)	45	29



	EFF
PSf	^{241}Am (%)
D (65 μm)	80





High EFF alpha emitters \rightarrow \varnothing 50 - 100 μm

Direct radionuclide determinations - PSm

Alpha emitter radionuclides.

^{241}Am



Direct radionuclide determinations - PSm


Beta emitter radionuclides.



<u>Particle size</u> <u>(μm)</u>	<u>^3H</u> <u>(%)</u>	<u>^{14}C</u> <u>(%)</u>	<u>$^{90}\text{Sr}/^{90}\text{Y}$</u> <u>(%)</u>
17	4,51	74,3	192,9
32	2,69	68,7	189,8
57	1,73	62,7	190,3
79	1,16	54,5	186,7
145	0,81	45,1	184,8

High EFF beta emitters \rightarrow ϕ 50- 100 μm

Low EFF low beta emitters ^3H




Direct radionuclide determinations - Calibration

Chemical Quenching.


Nitromethane : 0 - 7.57 $\mu\text{l}/\text{mL}$



^3H PSm 400 – 500 μm



^{36}Cl PSm 120 – 180 μm




EFF \downarrow \rightarrow Quenching parameter (SQP) \downarrow


Direct radionuclide determinations - Calibration

Colour Quenching.

Methyl Orange : 0 – 0.014 g/L



^{36}Cl PSm 120 – 180 μm



EFF \downarrow \rightarrow Quenching parameter (SQP) \downarrow



Direct radionuclide determinations - Calibration

Particle Quenching.

BaCl₂ : 0 – 1.4 M



³H PSm 120 – 180 μm



SQP 671 ± 3



Particle Quenching ↑ EFF ↓
Spectra position ~ slightly change
EFF ↓ → Quenching Parameter (density) ↑

Direct radionuclide determinations - Calibration


Optical Quenching.




Diameter (μm)	^{36}Cl (%)	SQP(E)
120-180	88.2(10)	680(9)
400-500	83.3(5)	737(4)




Low particle quenching



High particle quenching



High optical quenching



Low optical quenching

EFF \downarrow \rightarrow Quenching parameter (SQP) \uparrow

Direct radionuclide determinations - Calibration

PS Quenching.

Chemical Quenching → Quenching parameter (SQP)

Colour Quenching → Quenching parameter (SQP)

Particle Quenching → Quenching parameter (density)

Optical Quenching → Constant for a defined measurement conditions

Direct radionuclide determinations - Calibration

Ciemat- Nist procedure.



Energy distribution spectrum $S(E')$ that reach PSm: Simulation of particle Quenching by Montecarlo.


Direct Radionuclide determ - Alpha Beta discrimination

Pulse Shape distribution.



Time pulse shape distribution of alpha emitters is longer than beta emitters.

LS alpha \gg beta -- PS alpha $>$ beta.





Characteristic Light Pulse Shapes of Alpha and Beta Pulses in LSC

- Delay depends on 3T concentration.
- PS energy spreads faster by polymeric chains.
- Secondary solvent (Nafthalene, Disopropylnaphthalene) delays signal


Direct Radionuclide determ - Alpha Beta discrimination

Pulse Shape distribution.

Quantulus detector



Beta Classification LS = PS
 Alpha Classification PS ↓ LS
 Alpha misclassified even PSA=0 (fast)



Direct Radionuclide determ - Alpha Beta discrimination


Pulse Shape distribution.

Quantulus detector




PSm: PPO+POPOP+ Naphtalene (120-150 μm)

$^{90}\text{Sr}/^{90}\text{Y}$



^{241}Am



Alpha Pulses delayed
Alpha Classification PS improve

Direct Radionuclide determ - Alpha Beta discrimination

Pulse Shape distribution.

Triathler detector


PSm(120-150 μm)



$^{90}\text{Sr}/^{90}\text{Y}$

^{241}Am

PPO+POPOP+ Naphtalene



Alpha / Beta discrimination


Direct Radionuclide determ - Alpha Beta discrimination

Pulse Shape distribution.

Triathler detector



PSm(120-150 μm) PPO+POPOP+ 2g Naphtalene



Integrating areas: Beta emitter: misclassification < 2% - EFF 90%
Alpha emitter: misclassification 1% - EFF 25%

Direct Radionuclide determinations - Applications

Applications



- Limited use for routine radionuclide determinations.
- Useful for specific determinations.

Direct Radionuclide determinations - Applications

Applications



High Salty matrices.

LS – phase separation

PS - stable


Determination of Radiotracers in oil reservoirs

Institute for Energy Technology (Norway)

Radiotracer: $S^{14}CN^-$

Matrix: $NaClO_4$ conc.

Calibration: colour and particle quenching



	AA (dpm/g)	MA (dpm/g)
LS – IFE	4.23	0.55
PS - UB	$4.32 \pm .016$	0.59 ± 0.04

Direct Radionuclide determinations - Applications

Applications

Continuous detection.

LS – measurement cell: sample + cocktail mixture

- unstable and reagents consumption, waste generation.

PS – measurement cell : sample + solid scintillator

- stable – EFF detection

. Chromatographic determinations

. Environmental determination (^3H)


Multiple sample analysis.

. Microplates (labelled samples)


Selective Radionuclide determinations - PSresin

Plastic Scintillating resins – Extractive Scintillating resins

- . Alpha and Beta spectra distribution + Scintillation →
no selective technique / spectra overlapping.
- . PS → solid platform to implement separation procedures



PSm




Plastic Scintillating resin (PSresin)


Selective Radionuclide determinations - PSresin

Plastic Scintillating resins – Extractive Scintillating resins

- . Alpha and Beta spectra distribution + Scintillation →
no selective technique / spectra overlapping.
- . PS → solid platform to implement separation procedures



PSm




Plastic Scintillating resin (PSresin)


Selective Radionuclide determinations - PSresin

Plastic Scintillating resins – Extractive Scintillating resins.

. **PSresin**: separation step + measurement preparation



PSresin



Selective Radionuclide determinations - PSresin

Plastic Scintillating resins – Extractive Scintillating resins.

. **PSresin**: separation step + measurement preparation

- . Time reduction
- . Reagents and man power reduction
- . Waste reduction
- . Development in progress
- . Cost ?



Selective Radionuclide determinations - PSresin

Preparation approaches.

. Approaches to incorporate the selective capability:


- **Immobilization**: selective extractant (solvent) coating the PSm or CPS.
- **Imprinted polymers**: selective cavity on the PSm or CPS
- **Covalent bounding**: selective extractant bounded on the PSm or CPS surface.

-

Selective Radionuclide determ - Applications

Immobilization.

^{90}Sr in Water



Extractant: 4,4'(5')-di-t butylcyclohexane 18-crown-6 1M Octanol

Separation conditions: LiNO_3 6 M

	Act (dpm)	Act calc (dpm)	Error (%)
Drinking water	8.02	8.18	1.94
	7.77	7.66	-1.44
	7.66	7.54	-1.51
Sea Water	7.88	8.06	2.29
	8.00	7.80	-2.56
	7.70	8.01	4.09
River Water	7.86	7.84	-0.28
	7.75	7.60	-1.94
	7.77	7.72	-0.61


Selective Radionuclide determ - Applications

Immobilization.

Radiotracer $S^{14}CN$ in oil reservoir.

Extractant: Aliquat 336

Separation conditions: water



Sample	Conductivity (mS cm ⁻¹)	TOC (mg L ⁻¹)	Activity (Bq L ⁻¹)	
			PS resin	IFE
1	52.1 ± 0.1	13.6 ± 0.3	2.89 ± 0.14	2.98 ± 0.09
2	51.5 ± 0.2	18.4 ± 0.3	2.01 ± 0.04	2.21 ± 0.06
3	50.8 ± 0.2	22.4 ± 0.7	1.00 ± 0.09	1.02 ± 0.07
4	51.3 ± 0.1	39.2 ± 0.6	1.42 ± 0.02	1.36 ± 0.05
5	51.5 ± 0.2	18.4 ± 0.7	2.66 ± 0.08	2.76 ± 0.08


Selective Radionuclide determ - Applications

Immobilization.

^{99}Tc in water and urine.

Extractant: Aliquat 336

Separation conditions: HCl 0.1M



Sample	Activity (dpm L ⁻¹)	Activity Calc (dpm L ⁻¹)	Error (%)
Sea Water	24,3	23,0	-5,3
Sea Water	24,3	25,1	3,3
Sea Water	24,2	22,8	-6,2
Urine	0,43	0.44	2,4
Urine	0,46	0,42	-6.5

Selective Radionuclide determ - Applications

Immobilization.

^{210}Pb in water.

Extractant: 4,4'(5')-di-t butylcyclohexane 18-crown-6 1M Octanol


Separation conditions: HNO_3 2M

Sample	Activity (dpm/mL)	Activity Calc (dpm/mL)	Deviation (%)
Ebro river	10,1	10,8	-7,0
	10,1	11,0	-9,2
	10,1	9,8	3,6
Subterranean water	10,9	11,7	-7,0
	11,4	11,4	0,2
	11,4	11,8	-4,1
Congost river	11,0	11,4	-4,2
	10,4	10,7	-3,3
	11,4	11,3	0,7

Selective Radionuclide determinations

Imprinted.

Synthesis in presence of a template.



- Selectivity based on steric impediment
- No breakthrough volume
- Use of generic extractants

In progress: PSresin ^{63}Ni


Selective Radionuclide determinations

Imprinted.


Synthesis in absence of template.

Polymer Free Volume – Porosity

SEM



PSm





PSm-Porous

Selective Radionuclide determinations





Imprinted.

Synthesis in absence of template.

SEM



AFM



In progress: ^{222}Rn

Selective Radionuclide determinations

Covalent Bounding.

Selective extractant bounded on the PS surface.

- No breakthrough volume

Studies **in progress**:

- CPS - ^{233}U

- CPS – ^{63}Ni

Selective Radionuclide determ - Applications

Automation.

- Routine analysis improves by using automated systems.
- Inclusion of PSresin cartridges fits this approach

Summary and Future challenges.

Summary.

- Alpha and Beta emitters can be determined by PS (Low beta emitters)
 - Calibration by using Quenching Parameters
 - Alpha/ Beta discrimination is possible.
 - Selective determinations are possible
-
- Selective Plastic Scintillating Resins.

Summary and Future challenges.

Future challenges.

- Selective Plastic Scintillating Resins.
 - . Procedures of extractant incorporation
 - . Controlled porous materials
- New Plastic Scintillating formats.
 - . Foils
 - . Monolites
- Automation.
- New application fields.
 - . Routine control
 - . Medical
 - . Decommissioning
 - . Emergency situation

Plastic Scintillators and related analytical proposals

2nd International Workshop on Plastic Scintillation in Practice

Barcelona, Spring 2018

José F. García, Alex Tarancón, Héctor Bagán



UNIVERSITAT DE
BARCELONA



Thanks for your attention.