

APPLICATION OF LSC FOR ^3H , ^{14}C , ^{36}Cl , ^{41}Ca AND ^{63}Ni DETERMINATION IN VARIOUS MATRIX FROM NUCLEAR WASTE

FROM RESEARCH TO INDUSTRY

cea den

R. Brennetot, M. Giuliani, C. Colin, E. Laporte,
S. Guégan, J. P. Degros, A. Masset, P. Deloffre,
C. Mougel

Den – Service d'Etudes Analytiques et de Réactivité des Surfaces (SEARS)

CEA, Université Paris-Saclay

F-91191, Gif sur Yvette, France

rene.brennetot@cea.fr




Context of the study



Beta emitters isolation

^3H , ^{14}C in graphite and concrete
 ^{36}Cl , ^{63}Ni , ^{41}Ca




Matrices validation and / or proficiency test results



Conclusion and prospects

- The National Radioactive Waste Management Agency (**ANDRA**) is in charge of the long-term management of radioactive waste in France.
 - Controlling the different types of wastes is essential for the waste management.



Various types of Nuclear wastes



Radionuclide	Acceptance limit (Bq.g ⁻¹) for LLW
³ H	2 × 10 ⁵
⁶⁰ Co	1.3 × 10 ⁸
¹³⁷ Cs	3.3 × 10 ⁵
¹⁴ C	9.2 × 10 ⁴
³⁶ Cl	5
⁶³ Ni	3.2 × 10 ⁶
⁵⁵ Fe	6.1 × 10 ⁹
Σ alpha-emitters	3.7 × 10 ³

- ³H : Half life 12,3 years
- ¹⁴C : Half life 5700 years
- ³⁶Cl : Half life 3,02 10⁵ years
- ⁶³Ni : Half life 98,7 years
- ⁴¹Ca : Half life 1 10⁵ years

For surface disposal :

Maximum activity limit : depends of the RN





Beta emitters measurement : **Liquid Scintillation Counting**



Most widely used for isolated RN

RN need to be isolated first from the entire matrix



Pyrolysis method for
 ^3H and ^{14}C in all matrix


^{36}Cl in graphite

cea den COMBUSTION METHOD FOR $^3\text{H}/^{14}\text{C}$ IN GRAPHITE

For graphite :

T° rising from amb to 950°C
 in 20 min
 Then 1H at 950 °C
 Under O₂

COMBUSTION



^3H trapping : HCl 0.2M (2 traps)
 ^{14}C trapping : carbosorb (2 traps)

Sample boat empty at the end


➔ Extraction yield : 100%

cea den CARBOWASTE RESULTS FOR $^3\text{H}/^{14}\text{C}$ IN GRAPHITE

Carbowaste : 5 g of irradiated graphite segmented 2 + 2 (+ 1) with 2 determinations



Addition of ^3H activity obtained in the 2 first traps, more than 95% in trap 1

Addition of ^{14}C activity obtained in the 2 last traps, more than 95% in trap 3



^3H activity : 83,320 Bq/g \pm 5%

^{14}C activity : 21,310 Bq/g \pm 5%



Very good agreement


Methodology :

Use of labelled H-3 and C-14 concrete

optimization of temperature rising program for maximum tritium recovery by comparison with calculation from concrete preparation

Preparation of labelled concrete :

Use of H-3 glucose and C-14 glucose in H₂O



Cement added with water and mixed activity check deposited in vessel



Drying of concrete for at least 3 weeks



Weighing of concrete and calculation of H-3 and C-14 activity

concrete Grinding to obtain powder



To optimize temperature program, use of fine powder (< 250 μm obtained after sieving)


Different temperature program and gas flow tested :




**Only program 4
for recovery > 80%**

Different factors has been studied :

- particle sizes
- amounts



cea den MATRIX VALIDATION : CONCRETE



Sample dissolution
(acid digestion)

Ni, Co carrier

Nitroso 2
naphthol citrate


Liq-Liq extraction
for Co

Liq-Liq extraction
For Ni-DMG

ICP-AES
Ni


LSC
 ^{63}Ni

Quality check every month of the quenching curve
With certified activity sealed source




cea den BETA EMITTERS SEPARATION : ^{36}Cl


Quality check every month
With certified activity sealed source




^{36}Cl control



IC
Quantum yield



CARBOWASTE RESULTS FOR ³⁶Cl AND ⁶³Ni IN GRAPHITE



³⁶Cl activity : 27 Bq/g ± 10%


⁶³Ni activity : 113,460 Bq/g ± 10%

Very good agreement


Inter comparison on real matrix :

Concrete, metals, graphite, ...

→ too rare



cea den BETA EMITTERS SEPARATION : ^{41}Ca



Dissolution step :

→ Depends of the matrix


Concrete, metals : **acid digestion**


LSC measurement :

difficult for **low** contents
 due to poor efficiency of ^{41}Ca

→ Counting time high

If Sr high content





High importance for beta emitters measurement in nuclear waste


LSC is the reference method



^3H and ^{14}C extraction yield by pyrolysis has been validated for various matrices :

Concrete : 85% (^3H) 90% (^{14}C)

graphite : 100% for ^3H and ^{14}C



Associated rising temperature program

And gas flow values

^{63}Ni and ^{36}Cl radiochemistry efficiency proven with carbowaste results (graphite)

Need of inter comparison test on real samples (concrete, metals,.....)

^{41}Ca radiochemistry elaborated and quenching curve on TRICARB 2900 realized

Work in progress for ^3H measurement in another matrix : metals, resins, oils



Rene BRENNETOT

Den – Service d'Etudes Analytiques et de
Réactivité des Surfaces (SEARS)
CEA, Université Paris-Saclay
F-91191, Gif sur Yvette, France
rene.brennetot@cea.fr

Commissariat à l'énergie atomique et aux énergies alternatives
Centre de Saclay | 91191 Gif-sur-Yvette Cedex
T. +33 (0)1 XX XX XX XX | F. +33 (0)1 XX XX XX XX

Etablissement public à caractère industriel et commercial | RCS Paris B 775 685 019